



Barrett LRI Set

Titanium Instruments

4-622

Barrett LRI Diamond Knife



0.17mm blade thickness




- Blade preset at 550 micron setting / extension
- 8 facet vertical cutting blade
- Retractable long handle
- Length 125mm

The diamond knife has a unique double sided footplate for lateral support on the diamond blade, with an opening behind the footplate to give good visibility.

Barrett Limbal Markers




view on top of marker



- 12mm diameter part radial mark of length 6.0mm (ref: 9-741-1) or 8.0mm (ref: 9-741-2)
- Two 2.0mm side marks
- Instruments uniquely coloured for identification Gold (ref: 9-741-1) or Blue (ref: 9-741-2)
- Length 114mm


9-741-1

6.0mm arc length on 12.0mm diameter



9-741-2

8.0mm arc length on 12.0mm diameter



Duckworth & Kent Ltd.

Terence House
7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254 Fax: +44 (0)1462 896288
Web: www.duckworth-and-kent.com
Email: info@duckworth-and-kent.com

www.duckworth-and-kent.com

