



Duckworth & Kent has been a pioneering engineering company since it was formed in 1959. With over 40 years experience in manufacturing ophthalmic, neuro and ENT titanium surgical instruments, the company understands the needs of the surgeon and prides itself on customer satisfaction.

All Duckworth and Kent surgical instruments are manufactured by Duckworth and Kent, using titanium alloy as the raw material. Titanium is well tolerated by tissue and therefore well suited to the requirements for surgical instrumentation, offering advantages such as being non-corrosive, non-magnetic, relatively lightweight and yet extremely durable. The use of titanium is complemented by the advanced manufacturing methods using CNC equipment to achieve the most demanding tolerances to a consistent standard. The final result is a range of titanium neurosurgical instruments of exceptional quality with fine, delicate tips and teeth.

*Duckworth & Kent's
Guarantee:*

D&K instruments are guaranteed for life against defects in materials and workmanship.

Instrument Repair.

Full instrument repair service, contact D&K. Ensure to send your instrument repair to D&K via a traceable shipping method. Decontamination certification is required with each instrument.

*Custom Instrument
Design and Development.*

D&K produces custom instruments to accommodate new or modified surgical techniques, contact D&K.

Quality Standards.

Certified and Registered ISO9001, EN46001, ISO13485 and CE Mark Approved, which demonstrates the highest level of commitment to quality – The standard for excellence.

The schematic line drawings and photographs in this catalogue are fully protected by copyright. No part of this Catalogue may be reproduced in any form without prior written permission of Duckworth & Kent Ltd. Duckworth and Kent reserves the right to make changes at any time, without notice, in product specifications and availability.

Schematic line drawings of instruments are not scaled and tip drawings are shown enlarged.

Duckworth & Kent Ltd.

Terence House
7 Marquis Business Centre
Royston Road
Baldock, Herts, SG7 6XL, England

Telephone: +44 (0) 1462 893254
Fax: +44 (0) 1462 896288
Web: www.duckworth-and-kent.com
E-Mail: info@duckworth-and-kent.com

Duckworth & Kent USA Ltd.

10805 Sunset Office Drive, Suite 207
St Louis, MO 63127-1026 USA

Telephone: (314) 822-5400
(800) 756-7773
Fax: (314) 822-5426
Web: www.duckworth-and-kent.com
E-Mail: sales@dkusa.org





Table of Contents

Curette.....	1
Diamond Knives	
<i>Straight Handle</i>	3
<i>Bayonet Handle</i>	4
Dissectors	5
Elevators	6
Elevators, Knives & Blade Holders	7
Forceps	
<i>Straight – Flat & Round Handles</i>	8
<i>Straight Bi-Polar – Flat & Round Handles</i>	9
<i>Bayonet – Flat Handles</i>	9
<i>Bayonet – Round Handles</i>	11
<i>Bayonet Bi-Polar – Flat & Round Handles</i>	13
<i>Extended Squeeze Handle Forceps</i>	14
Needle Holders	
<i>Straight – Flat & Round Handles</i>	15
<i>Bayonet – Flat Handles</i>	16
<i>Bayonet – Round Handles</i>	17
Nerve Hooks & Probes	18
Probes & Spatulas	19
Scissors	
<i>Straight – Flat & Round Handles</i>	20
<i>Bayonet – Flat Handles</i>	21
<i>Bayonet – Round Handles</i>	24
Anastomosis	
<i>Scissors & Forceps</i>	25
<i>Needle Holders & Hooks</i>	26
Index	27


CURETTE

Angled Up Curette

NLS 5092



- Straight handle
- 2.25mm diameter, angle 45°
- Length 227.0mm




Angled Up Bayonet Curette

NLS 5092-1



- Bayonet handle
- 2.25mm diameter, angle 45°
- Length 224.0mm



Angled Up Bayonet Curette (short)

NLS 5098-1



- Bayonet handle
- 3.0mm diameter, angled 90°
- Length 255.0mm



Angled Down Bayonet Curette (short)

NLS 5098-2



- Straight handle
- 3.0mm diameter, angled 40°
- Length 255.0mm




Angled Up Bayonet Curette (short)

NLS 5098-3



- Straight handle
- 2.25mm diameter, angled 45°
- Length 255.0mm



CURETTE

Angled Up Bayonet Curette (long)

NLS 5099-1



- Bayonet long handle
- 3.0mm diameter angled 90°
- Length 280.0mm



Angled Down Bayonet Curette (long)

NLS 5099-2



- Bayonet long handle
- 3.0mm diameter angled 40°
- Length 280.0mm




Angled Up Bayonet Curette (long)

NLS 5099-3



- Bayonet long handle
- 2.25mm diameter angled 45°
- Length 280.0mm



Angled Up Bayonet Curette (long)

NLS 5099-4



- Bayonet long handle
- 3.0mm diameter angled 45°
- Length 280.0mm



Angled Up Bayonet Curette (long)

NLS 5099-5




- Bayonet long handle
- 3.0mm diameter angled 90°
- Length 280.0mm



DIAMOND KNIVES (Straight Handle)

Straight Neuro Diamond Knife (45° single edge blade)

NLS 4-080




- Retractable short handle
- Length 177.5mm



Straight Neuro Diamond Knife (45° double edge, vertical cutting blade)

NLS 4-090




- Retractable short handle
- Length 177.5mm



Straight Neuro Diamond Knife (45° single edge blade)

NLS 4100




- Retractable long handle
- Length 218.5mm



Straight Neuro Diamond Knife (45° double edge, vertical cutting blade)

NLS 4105




- Retractable long handle
- Length 218.5mm



Straight Neuro Diamond Knife (Tri-facet blade)

NLS 4115




- Retractable long handle
- Length 218.5mm




DIAMOND KNIVES (Bayonet Handle)

Bayonet Neuro Diamond Knife (45° single edge blade)

NLS 4400




- Bayonet handle with blade guard
- Length 207.5mm




Bayonet Neuro Diamond Knife (45° double edge, vertical cutting blade)

NLS 4405




- Bayonet handle with blade guard
- Length 207.5mm




Bayonet Neuro Diamond Knife (Tri-facet blade)

NLS 4415




- Bayonet handle with blade guard
- Length 207.5mm




Bayonet Neuro Diamond Knife (1.4mm lance blade)

NLS 4420




- Bayonet handle with blade guard
- Length 207.5mm




Bayonet Neuro Diamond Knife (30° single edge blade)

NLS 4425




- Bayonet handle with blade guard
- Length 207.5mm




Bayonet Neuro Diamond Knife (1.4mm spear blade)

NLS 4430




- Bayonet handle with blade guard
- Length 207.5mm



DISSECTORS

Bayonet Dissector (small)


NLS 5080



- Bayonet handle
- Length 230.0mm
- 1.15mm wide

Bayonet Dissector (medium)


NLS 5082



- Bayonet handle
- Length 230.0mm
- 1.7mm wide

Bayonet Dissector (large)

NLS 5084




- Bayonet handle
- Length 230.0mm
- 2.1mm wide

ELEVATORS

Elevator


NLS 5094



- Straight handle
- Length 227.0mm
- 2.0mm wide

Bayonet Elevator


NLS 5094-1



- Bayonet handle
- Length 224.5mm
- 2.0mm wide

Elevator


NLS 5095



- Straight handle
- Length 227.0mm
- 2.6mm wide

Bayonet Elevator

NLS 5095-1




- Bayonet handle
- Length 224.5mm
- 2.7mm wide


ELEVATORS, KNIVES & BLADE HOLDERS

Elevator (oval tip)

NLS 5097



- Straight handle
- Length 227.0mm
- 2.75mm wide




VIEW ON ARROW



SECTION A-A

Bayonet Elevator (oval tip)

NLS 5097-1



- Bayonet handle
- Length 224.5mm
- 2.75mm wide




VIEW ON ARROW



SECTION A-A

Bayonet Knife


NLS 5098-8



- Bayonet handle
- Length 255.0mm

Bayonet Blade Holder

NLS 5098-9




- Bayonet handle
- Length 255.0mm

VIEW ON ARROW


Straight Flat Handle Forceps

Straight Neuro Forceps (blunt tip)

NLS 4000




- 1.0mm Tip
- Straight flat handle
- Length 155.0mm




Straight Neuro Forceps (fine tip)

NLS 4002




- 0.5mm Tip
- Straight flat handle
- Length 155.0mm




Straight Round Handle Forceps

Straight Neuro Forceps (blunt tip)

NLS 5000




- 1.0mm Tip
- Straight round handle
- Length 180.0mm




Straight Neuro Forceps (fine tip)

NLS 5002




- 0.5mm Tip
- Straight round handle
- Length 180.0mm




Straight Flat Handle Forceps

Bi-Polar Straight Neuro Forceps (blunt tip)

NLS 4007




- Bi-Polar 1.0mm Tip
- Straight flat handle
- Length 155.0mm




Straight Round Handle Forceps

Bi-Polar Straight Neuro Forceps (blunt tip)

NLS 5007




- Bi-Polar 1.0mm Tip
- Straight round handle
- Length 190.0mm




Bayonet Flat Handle Forceps

Bayonet Neuro Forceps (straight, blunt tip)

NLS 4010




- Straight 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 230.0mm




Short Bayonet Neuro Forceps (straight, blunt tip)

NLS 4010-1




- Straight 1.0mm Tip
- Working length 63.5mm (2½ inch)
- Bayonet flat handle
- Length 208.0mm




Bayonet Neuro Forceps (straight, fine tip)

NLS 4011



- Straight 0.5mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 230.0mm




NLS 4011-1 *Short Bayonet Neuro Forceps (fine tip)*


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (straight, fine tip)

NLS 4012



- Straight 0.5mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 230.0mm




NLS 4012-1 *Short Bayonet Neuro Forceps (straight fine tip)*


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (curved up, blunt tip)

NLS 4014



- Curved up 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 230.0mm




NLS 4014-1 *Short Bayonet Neuro Forceps (curved up, blunt tip)*


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (curved up, fine tip)

NLS 4016



- Curved up 0.5mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 230.0mm




NLS 4016-1 *Short Bayonet Neuro Forceps (curved up, fine tip)*

- Working length 63.5mm (2½ inch)
- Length 208.0mm


FORCEPS (Bayonet Flat Handle & Bayonet Round Handle)

Bayonet Neuro Forceps (angled up, blunt tip)

NLS 4018



- Angled up 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 230.0mm




NLS 4018-1 *Short Bayonet Neuro Forceps (angled up, blunt tip)*


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (angled up, fine tip)

NLS 4020



- Angled up 0.5mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 230.0mm




NLS 4020-1 *Short Bayonet Neuro Forceps (angled up, fine tip)*

- Working length 63.5mm (2½ inch)
- Length 208.0mm


Bayonet Round Handle Forceps

Bayonet Neuro Forceps (straight, blunt tip)

NLS 5010



- Straight 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet round handle
- Length 230.0mm




NLS 5010-1 *Short Bayonet Neuro Forceps (straight, blunt tip)*


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (straight, fine tip)

NLS 5012



- Straight 0.5mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet round handle
- Length 230.0mm




NLS 5012-1 *Short Bayonet Neuro Forceps (straight, fine tip)*


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (curved up, blunt tip)

NLS 5014



- Curved up 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet round handle
- Length 230.0mm




NLS 5014-1 Short Bayonet Neuro Forceps (curved up, blunt tip)


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (curved up, fine tip)

NLS 5016



- Curved up 0.5mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet round handle
- Length 230.0mm




NLS 5016-1 Short Bayonet Neuro Forceps (curved up, fine tip)


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (angled up, blunt tip)

NLS 5018



- Angled up 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet round handle
- Length 230.0mm




NLS 5018-1 Short Bayonet Neuro Forceps (angled up, blunt tip)


- Working length 63.5mm (2½ inch)
- Length 208.0mm

Bayonet Neuro Forceps (angled up, fine tip)

NLS 5020




- Angled up 0.5mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet round handle
- Length 230.0mm



NLS 5020-1 *Short Bayonet Neuro Forceps (angled up, fine tip)*


- Working length 63.5mm (2½ inch)
- Length 208.0mm




Bayonet Flat Handle Bi-Polar Forceps

Bi-Polar Bayonet Neuro Forceps (straight, blunt tip)

NLS 4037




- Straight Bi-Polar 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet flat handle
- Length 220.0mm




Bayonet Round Handle Bi-Polar Forceps

Bi-Polar Bayonet Neuro Forceps (straight, blunt tip)

NLS 5037



- Straight Bi-Polar 1.0mm Tip
- Working length 89.0mm (3½ inch)
- Bayonet round handle
- Length 220.0mm




FORCEPS (Extended Squeeze Handle)

Extended Squeeze Handle Forceps (fine tip)

NLS 5030



- Squeeze handle
- Length 155.0mm




Extended Squeeze Handle Forceps (blunt tip)

NLS 5030-1




- Squeeze handle
- Length 155.0mm




Extended Squeeze Handle Forceps (fine tip, toothed)

NLS 5030-2



- Squeeze handle
- Length 155.0mm




NEEDLE HOLDERS (Straight - Flat & Round Handles)


Straight Flat Handle Needle Holder

Straight Needle Holder (straight jaw)

NLS 4060



- Straight jaw
- Straight flat handle
- Length 175.0mm




Straight Needle Holder (curved jaw)

NLS 4062




- Curved jaw
- Straight flat handle
- Length 175.0mm




Straight Round Handle Needle Holder

Round Handle Straight Needle Holder (straight jaw)

NLS 5060




- Straight jaw
- Straight round handle
- Length 180.0mm




Round Handle Straight Needle Holder (curved jaw)

NLS 5062




- Curved jaw
- Straight round handle
- Length 180.0mm




Flat Handle Bayonet Needle Holders

Bayonet Needle Holder (fine, straight jaw)

NLS 4064




- Straight jaw
- Bayonet flat handle
- Length 230.0mm




Bayonet Needle Holder (straight jaw with lock)

NLS 4065




- Straight jaw
- Bayonet flat handle with lock
- Length 230.0mm




Bayonet Needle Holder (fine, curved jaw)

NLS 4066



- Curved jaw
- Bayonet flat handle
- Length 230.0mm




NEEDLE HOLDERS (Bayonet - Round Handles)


Bayonet Round Handle Needle Holders

Bayonet Needle Holder (fine, straight jaw)

NLS 5064




- Straight jaw
- Bayonet round handle
- Length 230.0mm




Bayonet Needle Holder (straight jaw with lock)

NLS 5065



- Straight jaw
- Bayonet round handle with lock
- Length 230.0mm




Bayonet Needle Holder (fine, curved jaw)

NLS 5066



- Curved jaw
- Bayonet round handle
- Length 230.0mm




NERVE HOOKS & PROBES

Nerve Hook (bulbous tip)

NLS 5086


- Bayonet handle
- Length 230.0mm



Nerve Hook

NLS 5088

- Bayonet handle
- Length 230.0mm
- 90° 0.5mm tip




Mirror Probe

NLS 5090




- Straight handle
- Length 226.0mm
- 5.0mm diameter




Ball Tipped Probe

NLS 5093




- Straight handle
- Length 226.5mm
- 0.8mm diameter




Ball Tipped Bayonet Probe

NLS 5093-1




- Bayonet handle
- Length 224.5mm
- 0.8mm diameter




PROBES & SPATULAS

Turned Up Probe

NLS 5091





- Straight handle
- Length 224.0mm
- 0.8mm diameter



Teardrop Bayonet Probe

NLS 5099-6



- Bayonet handle
- Length 280.0mm



Spatula

NLS 5083



- Bayonet handle
- Length 225.0mm
- 2.0mm wide



Bayonet Spatula

NLS 5096-1



- Bayonet handle
- Length 224.5mm
- 2.4mm wide



Spatula

NLS 5099-7


- Bayonet handle
- Length 280.0mm
- 90° Tip




Straight Flat Handle Scissors

Straight Scissors (fine, straight blade)

NLS 4040




- Serrated upper blade
- Fine straight, 11.0mm long blades
- Straight flat handle
- Length 175.0mm




Straight Scissors (fine, curved blade)

NLS 4042




- Serrated upper blade
- Fine curved, 11.0mm long blades
- Straight flat handle
- Length 175.0mm




Straight Round Handle Scissors

Straight Scissors (fine, straight blade)

NLS 5040




- Serrated upper blade
- Fine straight, 11.0mm long blades
- Straight round handle
- Length 180.0mm




Straight Scissors (fine, curved blade)

NLS 5042




- Serrated upper blade
- Fine curved, 11.0mm long blades
- Straight round handle
- Length 180.0mm




Flat Handle Bayonet Scissors

Bayonet Scissors (regular, straight blade)

NLS 4044




- Serrated upper blade
- Fine straight, 11.0mm long blades
- Blunt tip
- Bayonet flat handle
- Length 230.0mm




Bayonet Scissors (regular, curved up blade)

NLS 4046




- Serrated upper blade
- Fine curved up, 11.0mm long blades
- Blunt tip
- Bayonet flat handle
- Length 230.0mm




Bayonet Scissors (fine, straight blade)

NLS 4048




- Serrated upper blade
- Fine straight, 11.0mm long blades
- Sharp tip
- Bayonet flat handle
- Length 230.0mm




Bayonet Scissors (fine, straight blade)

NLS 4049




- Serrated upper blade
- Fine straight, 11.0mm long blades
- Sharp tip
- Bayonet flat handle
- Length 218.0mm




Bayonet Scissors (fine, curved up blade)

NLS 4050




- Serrated upper blade
- Fine curved up, 11.0mm long blades
- Sharp tip
- Bayonet flat handle
- Length 230.0mm




Bayonet Scissors (fine, curved blade)

NLS 4051




- Serrated upper blade
- Fine curved up, 11.0mm long blades
- Sharp tip
- Bayonet flat handle
- Length 218.0mm




SCISSORS (Bayonet - Flat Handles cont.)

Bayonet Scissors (fine, straight blade)

NLS 4053




- Serrated upper blade
- Fine straight, 14.0mm long blades
- Sharp tip
- Bayonet flat handle
- Length 230.0mm




Bayonet Scissors (fine, curved blade)

NLS 4055




- Serrated upper blade
- Fine straight, 14.0mm long blades
- Sharp tip
- Bayonet flat handle
- Length 230.0mm




Bayonet Scissors (large blade)

NLS 4058




- Serrated upper blade
- Fine straight, large blades
- Pivot to tip length 24.5mm
- Blunt tip
- Bayonet flat handle
- Length 225.0mm




Round Handle Bayonet Scissors

Bayonet Scissors (regular, straight blade)

NLS 5044




- Serrated upper blade
- Fine straight, 11.0mm long blades
- Blunt tip
- Bayonet round handle
- Length 230.0mm




Bayonet Scissors (regular, curved up blade)

NLS 5046




- Serrated upper blade
- Fine curved up, 11.0mm long blades
- Blunt tip
- Bayonet round handle
- Length 230.0mm




Bayonet Scissors (fine, straight blade)

NLS 5048




- Serrated upper blade
- Fine straight, 11.0mm long blades
- Sharp tip
- Bayonet round handle
- Length 230.0mm




Bayonet Scissors (fine, curved up blade)

NLS 5050




- Serrated upper blade
- Fine curved up, 11.0mm long blades
- Sharp tip
- Bayonet round handle
- Length 230.0mm




Round Handled Scissors

Round Handle Scissors (fine, curved)

AN 1010




- Fine curved, 10.0mm long blades
- Round handle, length 160.0mm




Round Handle Scissors (curved)

AN 1011




- Curved, 18.0mm long blades
- Round handle, length 160.0mm




Round Handled Forceps

Straight Dissecting Round Handle Forceps (fine tips)

AN 1020




- Fine tips
- Round handle, length 160.0mm




Curved Dissecting Round Handle Forceps (fine tips)

AN 1021




- Fine tips
- Round handle, length 160.0mm




Straight Dissecting Round Handle Forceps (blunt tips)

AN 1022




- Blunt tips
- Round handle, length 160.0mm




Angled Dissecting Round Handle Forceps (fine tips)

AN 1023




- Blunt tips
- Round handle, length 160.0mm




Round Handled Needle Holders

Round Handle Needle Holder (fine, curved jaw)

AN 1030




- Fine curved jaw
- Round handle, length 160.0mm




Round Handle Needle Holder (curved jaw)

AN 1031




- Curved jaw
- Round handle, length 160.0mm



Hooks

Hook (90° angled tip)

AN 1060




- 90° angle tip
- Round handle, length 160.0mm



Hook (90° angled, bulbous tip)

AN 1061




- 90° bulbous tip
- Round handle, length 160.0mm



Hook (hooked tip)

AN 1062



- Hooked tip
- Round handle, length 160.0mm





INDEX

AN 1010.....	25	NLS 4018-1.....	11	NLS 4400.....	4	NLS 5040.....	20	NLS 5095.....	6
AN 1011.....	25	NLS 4020.....	11	NLS 4405.....	4	NLS 5042.....	20	NLS 5095-1.....	6
AN 1020.....	25	NLS 4020-1.....	11	NLS 4415.....	4	NLS 5044.....	24	NLS 5096-1.....	19
AN 1021.....	25	NLS 4037.....	13	NLS 4420.....	4	NLS 5046.....	24	NLS 5097.....	7
AN 1022.....	25	NLS 4040.....	20	NLS 4425.....	4	NLS 5048.....	24	NLS 5097-1.....	7
AN 1023.....	25	NLS 4042.....	20	NLS 4430.....	4	NLS 5050.....	24	NLS 5098-1.....	1
AN 1030.....	26	NLS 4044.....	21	NLS 5000.....	8	NLS 5060.....	15	NLS 5098-2.....	1
AN 1031.....	26	NLS 4046.....	21	NLS 5002.....	8	NLS 5062.....	15	NLS 5098-3.....	1
AN 1060.....	26	NLS 4048.....	21	NLS 5007.....	9	NLS 5064.....	17	NLS 5098-8.....	7
AN 1061.....	26	NLS 4049.....	22	NLS 5010.....	11	NLS 5065.....	17	NLS 5098-9.....	7
AN 1062.....	26	NLS 4050.....	22	NLS 5010-1.....	11	NLS 5066.....	17	NLS 5099-1.....	2
NLS 4000.....	8	NLS 4051.....	22	NLS 5012.....	11	NLS 5080.....	5	NLS 5099-2.....	2
NLS 4002.....	8	NLS 4053.....	23	NLS 5012-1.....	11	NLS 5082.....	5	NLS 5099-3.....	2
NLS 4007.....	9	NLS 4055.....	23	NLS 5014.....	12	NLS 5083.....	19	NLS 5099-4.....	2
NLS 4010.....	9	NLS 4058.....	23	NLS 5014-1.....	12	NLS 5084.....	5	NLS 5099-5.....	2
NLS 4010-1.....	9	NLS 4060.....	15	NLS 5016.....	12	NLS 5086.....	18	NLS 5099-6.....	19
NLS 4011.....	10	NLS 4062.....	15	NLS 5016-1.....	12	NLS 5088.....	18	NLS 5099-7.....	19
NLS 4011-1.....	10	NLS 4064.....	16	NLS 5018.....	12	NLS 5090.....	18		
NLS 4012.....	10	NLS 4065.....	16	NLS 5018-1.....	12	NLS 5091.....	19		
NLS 4012-1.....	10	NLS 4066.....	16	NLS 5020.....	13	NLS 5092.....	1		
NLS 4014.....	10	NLS 4-080.....	3	NLS 5020-1.....	13	NLS 5092-1.....	1		
NLS 4014-1.....	10	NLS 4-090.....	3	NLS 5030.....	14	NLS 5093.....	18		
NLS 4016.....	10	NLS 4100.....	3	NLS 5030-1.....	14	NLS 5093-1.....	18		
NLS 4016-1.....	10	NLS 4105.....	3	NLS 5030-2.....	14	NLS 5094.....	6		
NLS 4018.....	11	NLS 4115.....	3	NLS 5037.....	13	NLS 5094-1.....	6		



NOTES

